



WARD CARE SERVICES

Excellent quality



Bed frame materials: four pieces of
-grade PP guardrails are used
rest adjustment: $45^{\circ} \pm 5^{\circ}$
-rest system: luxury centrally-con-
ed casters for one-foot brake



Five-function electric bed

Model: CH-DH-III

PRODUCT FEATURES

- ◇ Bed specifications: L2150 x 1080 x 460/680m
- ◇ Backrest adjustment: $80^{\circ} \pm 5^{\circ}$
- ◇ Forward and backward tilt angles: $8^{\circ} \pm 1^{\circ}$
- ◇ luxury centrally-controlled brake, one brake for overall stopping
- ◇ Load-bearing capacity of bed $\geq 260\text{kg}$

- ◇ Overall lifting: 460-680mm
- ◇ Legrest adjustment: 45°±5°
- ◇ Four European-style casters that can be turned up and down
- ◇ The bed board is formed by cold-rolled steel with slip resistance and air

- ◇ Overall lifting: 460-680mm
- ◇ Legrest adjustment: 40° ±5°
- ◇ luxury centrally-controlled brake, one brake for overall stopping
- ◇ Load-bearing capacity of bed ≥260kg



CH-A2

Double-crank care bed

Bed specification: L2150x1080x500mm



CH-B1

Three crank care bed

Bed specifications: L2150 x 900 x 430/630mm



CH-B2

Double-crank care bed

Bed specifications: L2150 x 950 x 500mm



CH-B3

Double-crank care bed (ABS cranks)

Bed specifications: L2150x 900 x 500mm



CH-C2

Double-crank care bed

Bed specifications: L2150 x 900 x 500mm





CH-C4

Punched double-crank bed with hanging ABS bed head

Bed specifications: L2150 x 900 x 500mm



hanging

500mm

CH-C6

Punched single-crank bed with hanging ABS bed head

Bed specifications: L2150 x 900 x 500mm




hanging

500mm



[CH-C8] Stripped crank bed with stainless-steel bed head

Product specifications:
L2020 x 900 x 500mm




[CH-C11] Stripped single-crank bed with steel bed head

Product specifications:
L2020 x 900 x 500mm



[CH-D3] Striped flat bed with stain-resistant-steel bed head

Product specifications:
L2020 x 900 x 500mm



[CH-C9] Stripped double-crank bed with steel bed head

Product specifications:
L2020 x 900 x 500mm



[CH-D1] Punched single-crank bed with hanging ABS bed head

Product specifications:
L2150 x 900 x 500mm



[CH-D4] Striped flat bed with steel bed head

Product specifications:
L2020 x 900 x 500mm



Stainless-steel infantette
Product specifications:
L1840 x W900mm




[A21] Stainless-steel tiltable infantette
Product specifications:
L1840 x 530 x 920mm



Colorful stainless-steel tiltable
Product specifications:
L1810 x 520 x 920mm



[A24] ABS-1 Baby carriage
Product specifications:
L810 x 520 x 910mm



Full-covered child care bed
Product specification: L1880 x W900 x H60mm



[CH-T3] Child ward bed
Product specification: L1950/2100 x 900mm /1030 x 500mm

Turning-over care bad

(Full bending) (optional wood bed head)



Model: CH-F1

PRODUCT FEATURES

- ◇ Bed specification: L1980 x 900 x 500mm (inner diameter)/L2110 x 500 x 900mm (outer diameter)
- ◇ Turning-over adjustment: 0-65°
- ◇ Back-lifting adjustment: 0-85°
- ◇ legrest adjustment: downward leg bending 0-65° upward leg bending 0-30°
- ◇ Back lifting board size: 200x270mm



Punched double-crank bed with composite board bed head



PRODUCT FEATURES

- ◇ Bed specifications: L2070 x W900 x H500mm
- ◇ Backrest adjustment: $0-80 \pm 5^\circ$
- ◇ Legrest adjustment: $0-40 \pm 5^\circ$


Punched double-crank bed with composite board bed head and fence



PRODUCT FEATURES

- ◇ Bed specifications: L2070 x W950 x H500mm
- ◇ Backrest adjustment: $0-80 \pm 5^\circ$
- ◇ Legrest adjustment: $0-40 \pm 5^\circ$





Luxurious wheeled stretcher (trolley)



Model: B2

PRODUCT FEATURES

- ◇ Bed specifications: L 1900x720 x500/800mm
- ◇ Lifting distance: 300mm
- ◇ Surface material: die-cast formed with ABS engineering injection molding materials with high fire resistance and corrosion resistance
- ◇ Back-lifting adjustment: 65°
- ◇ Lifting adjustment: 500-800mm
- ◇ Brake system: 125mm luxury silent casters with flexible rotation and strong stability



[CH-Q1] Stainless-steel mulberry
garden-style gynecological bed



[CH-Q2] Stainless-steel porta-
-ble obstetrics table




[B4] Stainless-steel stretcher with two
big and two small wheels (size
customizable)

Product specifications:
L 1950 x 600 x 800mm



Ambulance folding stretcher

Highest product specification: L 1900
550 x 920mm
Lowest product specification: L 1900
550 x 250mm



[B5] Stainless-steel treatment trolley
(size customizable)

Product specifications:
L 650 x 450 x 800mm



[B5] Stainless-steel three-tier treatment trolley
(size customizable)

Product specifications:
L 750 x 450 x 800mm



[B10] **Stainless-steel fan-shaped instrument trolley** (size customizable)
Product specifications:
L600/450 x 800mm



[B10] **Stainless-steel three-tier instrument trolley** (size customizable)
Product specification: large: L750x 450 x 900mm
/L670x 450 x 900mm
/L600x 420 x 900mm



[B13] **Stainless-steel emergency trolley** (size customizable)
Product specifications:
L600 x 850mm



[B13] **Stainless-steel fully-enclosed anesthetic trolley** (size customizable)
Product specifications:
L650x 450 x 850mm




[B16] **Stainless-steel medication trolley** (size customizable)
Product specifications:
L700 x 900mm



[B16] **Stainless-steel infusion trolley**
Product specifications:
L700 x 450 x 900mm



[B17] **Stainless-steel enclosed delivery trolley** (size customizable)
Product specifications:
L1000 x 550 x 900mm



[B20] **Stainless-steel debris trolley** (size customizable)
Product specifications:
L800 x 350 x 700mm



[B18] **Stainless-steel delivery trolley** (size customizable)
Product specifications:
L1000 x 550 x 900mm



[B21] **Stainless-steel single-rod square-plate tray** (size customizable)
Product specifications:
L500 x 400 x 800 / 1200mm



[B23] **Stainless-steel laundry bag trolley** (size customizable)
Product specifications:
L900 x 700 x 950mm



[B25] **Stainless-steel care trolley** (size customizable)
Product specifications:
L950 x 550 x 800mm